

SURE CHECK® HIV Self-Test



SURE CHECK® HIV Self-Test

SURE CHECK® HIV Self-Test

One Drop of Blood is All You Need



A single-use in-vitro diagnostic self-test that detects antibodies to HIV-1 and HIV-2 in fingerstick whole blood. The test is intended for use by untrained lay users in a private setting as a self-test to aid in the diagnosis of HIV infection.



EASY TO USE

- Quick blood collection; needs only 2.5 µl of blood
- Unique Barrel Design: Allows direct transfer of blood from finger to test device



CONVENIENT

- Low hands on time, no bottles, no mixing, no pouring
- Innovative test design prevents accidental spilling



TRUSTED & PROVEN

- Over 13 million tests used professionally
- High Accuracy: Sensitivity: 97%; Specificity: 100%*

*Refer to product insert for more details

MOST COMPLETE KIT COMPONENTS COMPARED TO LEADING BRANDS

CONTENTS IN THE BOX



CONTENTS IN THE POUCH



SUITABLE FOR

- Individuals who believe they may have had exposure to the HIV virus or a risk event within the last 3 months period of time who are ages 14 and older.

RISK EVENTS

- Individuals who have sex with someone who is HIV positive or whose status you do not know
- Sexually active individuals, regardless of their sexual orientation or gender
- Individuals who use shared needles or syringes
- Individuals who perceive themselves as at risk and have never undergone HIV screening or been diagnosed with HIV

TEST FOR HIV IN 3 EASY STEPS WITH SURE CHECK® HIV Self-Test

1 PREPARE

FILL TIP OF TESTING DEVICE WITH BLOOD

2 RUN

POSITION TEST DEVICE VERTICALLY ABOVE TEST STAND

1st SNAP
2nd SNAP
3rd SNAP

PUSH HARD THROUGH THE FOIL YOU WILL FEEL 3 SNAPS

3

CHECK FOR PINK/PURPLE FLOW IN 1 MINUTE

PINK/PURPLE FLOW ✓
NO PINK/PURPLE FLOW ✗

IF NO PINK/PURPLE FLOW, PUSH TEST DEVICE FURTHER DOWN

4 READ RESULTS

WAIT 15 MINUTES

DO NOT READ THE RESULTS AFTER 20 MINUTES

NEGATIVE RESULTS

← CONTROL
TEST

C T

If test shows **ONLY** the Control Line, your test result is **NEGATIVE**.

Test again in 3 months.

POSITIVE RESULTS

FAINT LINE ← CONTROL → TEST ← DARK LINE

C T

If test shows **BOTH** the Control Line and Test Line, your test result is **POSITIVE**.

Visit your nearest HIV Testing Center or Health Care Professional.

INVALID RESULTS

← CONTROL → TEST

C T

If your results do not match negative or positive, the result is **INVALID**. It is necessary to run another test.

Visit your nearest HIV Testing Center or Health Care Professional.

SCAN ME! VIDEO

Ordering Information

Cat. No.	Product Description	Test Specimen*	Product Format	Intended Use	Intended Use	Storage Temp.	Packing Size
60-9508-M	SURE CHECK® HIV Self-Test	WB	Sampler with buffer vial	Detection of antibodies to HIV Type 1 and 2	For Self-Testing	8-30°C	1 T/ Kit

Remarks: *WB – Whole Blood

Distributed by

Reszon.
Diagnostics
INTERNATIONAL

Manufactured by

CHEMBIO
DIAGNOSTICS MALAYSIA

Reszon Diagnostics International Sdn. Bhd. 201001033978 (917901-W)

✉ info@reszonics.com @ www.reszonics.com

Doc No.: MM-02-03-08-EN
Rev. No.: 01
Effective Date: 2024.03.14

testnow✓



Reszon Diagnostics