



TEST YOURSELF,
PROTECT OTHERS

SALIXIUM COVID-19 RAPID ANTIGEN TEST (SALIVA/NASAL SWAB SAMPLES)

REF | SLXHB1-0621001 | ARTG Identifier: 418208



The non-intrusive COVID-19 rapid test kit for SELF-TEST in human saliva and nasal swab samples.

Get your results at ease within 15 minutes from the comfort of your home or any premise of your choice.



PRODUCT SPECIFICATION



Sensitivity*
97.24%



Test time
15 minutes



Sample type
Saliva and anterior
nasal swab specimen



Specificity*
100%



Storage condition
4-30°C



Shelf life
24 months

* vs FDA Emergency Use
Authorized real-time RT-PCR

SALIXIUM HIGHLIGHTS



Fast
Rapid screening in less
than 15 minutes



No Pain
Less invasive sample collection
method (saliva & nasal swab)



Halvec Verified
Tested by Halvec Laboratories,
free from alcohol and animal DNA



Safe Swabs
Disposable swabs for saliva and
nasal are included in individual
sterile packing to minimize
cross contamination




Biohazard Waste Bag
Disposable biohazard waste bag
is included for safe and easy
disposal after use

AUSTRALIAN SPONSOR / DISTRIBUTOR

Avantua Pharmaceuticals Pty Limited
7, 71 Queens Rd Melbourne 3004 Australia
P.O. BOX 6009 Melbourne 3004 Vic Australia

☎ 1300 418 188
@ www.avantua.com.au
✉ support@avantua.com.au



Doc. No.: MM-04-04-11-EN-AU
Rev. No.: 02
Effective Date: 2024.09.10